

Overview of Windows support

Product	Version	Windows XP (SP 3)	Windows Vista 32 bit (SP1)	Windows 7 32 bit (SP 1)	Windows 7 64 bit (SP 1)	Windows 8.0 (64 bit)	Windows 8.1 (64 bit)	Windows 10 (64 bit)
The Observer XT	10.0*, 10.1*	x		x				
	10.5*	x		x	x			
	11*			x	x			
	11.5*			x	x	x		
	12.0				x		x	
	12.5				x		x	
	13.0				x			x
	14.0					x		x
EthoVision XT	8*	x		x				
	8.5*	x		x	x			
	9*			x	x			
	10				x		x	
	10.1				x		x	
	11				x		x	
	11.5				x		x	x
	12				x			x
CatWalk	7*	x						
	8*	x	x					
	9.1	x		x	x			
	10			x	x			
	10.5			x	x			
	10.6				x		x	x
FaceReader	1*	x	x					
	2*	x	x					
	3*	x	x	x				
	4*	x		x	x			
	5*			x	x			
	5.1				x	x		
	6				x		x	
	6.1				x		x	x
7				x			x	
DanioScope	1				x		x	x
	1.1				x			x
uLog	2*	x						
	3.5				x			x

*This software is no longer supported.

Product	Version	Windows XP (SP 3)	Windows Vista 32 bit (SP1)	Windows 7 32 bit (SP 1)	Windows 7 64 bit (SP 1)	Windows 8.0 (64 bit)	Windows 8.1 (64 bit)	Windows 10 (64 bit)
MPEG Recorder	1*	x						
Theme	4	x						
	5	x						
ULTRAVOX	2*	x						
	2*	x						
ULTRAVOX XT	3				x			
	3.1				x		x	x
Media Recorder	1*	x		x				
	2*			x	x			
	2.5				x	x		
	3				x		x	
	4				x			x
Viso	1				x		x	
	2				x		x	
	3						x	
	4							x
ERASMUS LADDER	1						x	
	1.1						x	
	2							x
SITE LICENSE MANAGER	1				x			x
	1.1				x			x

*This software is no longer supported.

Details

We test and support the most recent service pack, with the US English professional or business version of Windows. Windows Vista is supported for some versions, but not recommended.

Other systems

We do not support Apple OS X (also not with emulators or similar), or Linux. Pocket Observer 3 runs on Windows Mobile 6.1 (version 3.1) and Android 2.2 and higher (Pocket Observer version 3.3).

Detailed specifications

For detailed system requirements of current Noldus software and hardware, please see the document, 'system requirements for Noldus software' on <http://www.noldus.com/downloads>.

Windows 10

For the products that are released after January 2016, we support Windows 10 (64 bit), unless specified otherwise. Note however, that for certain hardware devices, there may always not be Windows 10 drivers available from their manufacturers. For the products released before January 2016, we tested whether the most recent version works with Windows 10. See the table whether your version is supported with this operating system.

Legal notice.

These specifications are not exhaustive and are subject to change without notice and no rights can be derived from this document. All products and services supplied by Noldus are subject to our general terms and conditions.

This document was last updated May 2017.

Specifications subject to change without notice. No rights may be derived from this document.

© Noldus Information Technology 2017