



PuriCare Bedding Disposal Station

Protects operators from airborne particulates during bedding disposal. Large access opening allows for free range of motion. HEPA filter traps allergens; activated carbon and alumina exhaust filter adsorbs odors, contaminants and other gases.

◀ **Labconco**
www.labconco.com

Indigo RO Reverse Osmosis

Modular water purification system removes 99.99% of most water contaminants to protect animal health. Offers a complete graphical view of system from the user-friendly touch-screen display. Automatic flushing reduces time, labor and energy.

Edstrom ▶
www.edstrom.com



AllerGard NU-619 Animal Transfer Station

New lightweight model increases mobility in the lab. 10" of knee space allow the user to sit comfortably. Service-friendly microprocessor controls increase reliability. Energy Saver motor option lowers noise and energy costs.

◀ **NuAire, Inc.**
www.nuaire.com

DanioVision

High-throughput tracking system for zebrafish larvae in a multi-well plate. Provides a controlled environment with an innovative optical design for a distortion-free image. Allows detailed tracking and automation of experiments.

Noldus Information Technology ▶
www.noldus.com



To submit product releases for consideration,
please send artwork and descriptive copy to: info@labanimal.com

For more information about the Product Showcase feature,
please call Daina Tilton at **212.451.8567**



Zebrafish Housing Racks

Tubular stainless steel welded rack with open frame design to allow maximum visibility without sacrificing structural integrity. Onboard water delivery and drainage system is easily removable for cleaning and maintenance.

◀ **Aquaneering, Inc.**
www.aquaneering.com

Zebrafish Tanks

Durable, stackable self-cleaning tanks feature a subtle "V" formed into the floor, channeling solids to the rear of the tank. Removable baffle forces water along the bottom, carrying solids out of the overflow port.

Aquaneering, Inc. ▶
www.aquaneering.com



design, integrity, reliability.



the possibilities are endless.



Internationally respected for its products, years of design and manufacturing experience make Lomir your choice for quality innovative products.

INFUSION SYSTEMS

TETHERED AND AMBULATORY SYSTEMS FOR ALL LAB SPECIES.

ANIMAL JACKETS

ALL ANIMALS, ALL APPLICATIONS AND ALL THE ACCESSORIES.

JACKET & UNDERSHIRT ENSEMBLES
 FOR EXTERNAL TELEMETRY SYSTEMS.

PROTECTIVE COLLARS

RESTRAINTS

SNUGGLE FAMILY AND SLING RESTRAINTS


ANIMAL HANDLING & PROTECTIVE EQUIPMENT

ENVIRONMENTAL ENRICHMENT PRODUCTS

We design and manufacture products to meet your specifications.
 For more information contact us to discuss your requirements.

877.425.3604
www.lomir.com





Vevo LAZR Photoacoustic Imaging Platform

Real-time in vivo imaging of deep tissue. High sensitivity and specificity and inherent co-registration for oncology research (sentinel lymph node, oxygenation, tumor and molecular imaging) and contrast applications (cell specificity, biomarker studies).

VisualSonics
www.visualsonics.com


Mouse Rat Blood Pressure

Revolutionary non-invasive blood pressure testing for mice, rats and other laboratory animals. A compact, simple yet versatile turn-key system can test from 1 to 200 animals at a time, with independent control of each channel.

IITC Life Science Inc.
www.iitcinc.com



AUTOMATED BLOOD SAMPLER



Instech's second-generation sampler, the ABS2™, collects blood samples from laboratory animals automatically, reducing animal stress and improving research efficiency. It works with mice, rats and large animals alike, deposits blood either into refrigerated vials or onto dried blood spot paper, and features advanced sensors to ensure successful collection of pure blood precisely at your programmed time points.

INSTECH SOLOMON
www.instechlabs.com
Tel 800-443-4227



pH Natural 5.5 Powder-Free Examination Gloves

Help maintain a lower skin pH, less dryness and less skin irritation than standard gloves. Patent-pending technology aids in restoring natural skin pH, while moisturizing dry, chapped hands. Available in latex and soft nitrile.

◀ **Shentex**
www.shentex.com

Intensive Care Unit with Dome Cover

Allow a full range of respiratory therapeutic procedures, including administering oxygen, humidity and nebulization in combination with heat. Small units accommodate animals up to 20 pounds; medium/large units, up to 80 pounds.

ThermoCare ▶
www.thermocare.com



The Green Vivarium Foundation

Solutions

Resources

Recognition

Together, we can build a Greener Vivarium

Become a sponsor or member today!

www.greenvivarium.org



Braintree Scientific

Providing essential products to the life science industry for over 30 years.



Now offering *Wheaton Science Glass Products*



- ▶ Ampules
- ▶ Glass Bottles
- ▶ Staining Ware
- ▶ Vials
- ▶ Cell Culture
- ▶ Chemistry Glassware
- ▶ Chromatography
- ▶ Liquid Handling
- ▶ Lyophilization
- ▶ Plastic Bottles
- ▶ Caps & Closures

We offer over 500 products. If you don't see what you need, we will do our best to find it or develop it.



BRAINTREE SCIENTIFIC, INC.
www.braintreesci.com

Telephone: 781-348-0768 • Fax: 781-843-7932 • Email: info@BraintreeSci.com