

Overview of Windows support

Product	Version	Windows 2000	Windows XP	Windows Vista 32 bit	Windows 7 32 bit (32 or 64 bit PC)	Windows 7 64 Bit
The Observer	5	x	x			
	6		x			
	7 - 9		x	x		
	10.1		x	x	x	
	10.5		x		x	x
EthoVision	2 - 3	x	x			
	4 - 5		x			
	6		x	x		
	7.1		x		x	
	8		x		x	
	8.5		x		x	x
CatWalk	7		x			
	8		x	x		
	9.1		x		x	x
FaceReader	1		x			
	2		x	x		
	3		x	x	x	
	4		x		x	x
uLog	2		x			
Theme	4		x			
	5		x			
MPEG Recorder	1		x			
	2		x			
Media Recorder	1		x		x	
	2				x	x

Details

We test and support the most recent service pack, with the US English professional or business version of Windows. Windows Vista is supported for some versions, but not recommended.

Other systems

Noldus software only runs on Windows systems, not on Apple OS X (also not with emulators or similar), Android, or similar.

Detailed specifications

For detailed system requirements of current Noldus software and hardware, please see the document, 'system requirements for Noldus software' on <http://www.noldus.com/downloads>.